

Manual Exam Table with Four Storage Drawers

Model 5240

UMF Medical Manual Exam Tables with pneumatic-assist back have a 500lb weight capacity, and reversible side and front storage drawers to accommodate table placement in any room configuration. The table also features hidden stirrups, and an all-steel OSHA-compliant step with anti-slip guards. Ships assembled for easy installation.



Key features at a glance

- › 500lb capacity (226.8kg)
- › Removable, Seamless Upholstered Top with Prefix® Protective Finish
- › Galvanized steel base, heavy gauge, powder coated
- › Comfort Grip Back Release Handle
- › One-touch, easy glide stirrups with 4 lateral positions and friction lock length adjustment
- › Slide out leg section with removable upholstered leg rest
- › Removable 9"x9" plastic drip pan (22.9cm x 22.9cm)
- › Pull-out all steel non slip step 7" high (17.8cm)
- › 21" Paper roll holder (53.3cm)
- › 2 front drawers on full extension roller bearing glides 18.25"x15.5"x4" each (46.4cm x 39.4cm x 10.2cm)
- › 2 reversible side storage drawers on full extension roller bearing glides 18.25"x15.5"x4" each (46.4cm x 39.4cm x 10.2cm)
- › 5-year limited warranty

Specifications

- › Height: 32.75" (83.2cm)
- › Length: 58.5" – 74.5" (148.6cm – 189.23cm)
- › Width: 28" (71.12cm)
- › Stirrups: Infinitely adjustable
- › Net Weight: 241lb (109kg)
- › Shipping Weight: 320lb (145kg)
- › Shipping Dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

Accessories

- › 021 Chair-to-bed Elevated Leg Pad for Models 5240/5250 creates flat transfer surface. Fits snugly over leg pad, secured w/ strap
- › 251 Bierhoff Knee Crutches (Pair)
- › 350 Pelvic Tilt
- › 365 IV Pole (Patient Left)
- › 366 IV Pole (Patient Right)
- › 375 Safety Rail Package
- › 464 Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Left) Compatible with WA Green Series Light 48810, Field Installable for 5200 Series
- › 465 Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Right) Compatible with WA Green Series Light 48810, Field Installable for 5200 Series
- › 820 Removable Stainless Steel Drain Pan
- › 822 Removable Stainless Steel Drain Pan w/Spout
- › 849 220v Drawer Warmer (Factory Installed) (International)
- › 851 110v Drawer Warmer (Factory Installed)
- › 852 110v Outlet Package (Factory Installed)
- › 855 220v Outlet Package (Factory Installed) (International)
- › Fire Retardant Foam (Cal133) Available

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2020 United Metal Fabricators Inc. dba UMF Medical 06.19 All Rights Reserved.