

Treatment Table with Adjustable Backrest

Model 5570

With a 350-lb weight capacity, the Treatment Table with Adjustable Backrest provides a sturdy, durable surface for treatments. A five-year warranty offers peace of mind while the comfortable upholstery with PreFixx® Protective Finish makes the table easy to clean and disinfect.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



Key features at a glance

- › 350lb (226.8kg) weight capacity
- › Seamless upholstered top with Prefixx® Protective Finish
- › Manual, Adjustable Backrest
- › Available in all upholstery color options
- › All-welded steel construction
- › Leveling glides
- › Paper cutter band
- › Paper rod holder

Specifications

Height: 31.5" (80cm)
 Length: 70"- 81" with leg section (177.8cm -205.74cm)
 Width: 26.5" (67.3cm)
 Net Weight: 105lb (47.6kg)
 Shipping Weight: 158lb (71.7kg)
 Shipping Dimensions: 76"L x 30"W x 41"H (193cm x 76.2cm x 104.1cm)

Accessories

015 – Foot Section and pad
 826 – Painted Steel Storage Shelf