

H-Brace Treatment Table with Two Drawers

Model 5585

With a 350-lb weight capacity, the H-Brace Treatment Table with Two Drawers provides a sturdy, durable surface for treatments. A five-year warranty offers peace of mind while the comfortable upholstery with PreFixx® Protecti Finish makes the table easy to clean and disinfect.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



Key features at a glance

- › 350lb (158.8kg) weight capacity
- › Seamless upholstered top with Prefixx® Protective Finish
- › All-welded steel construction
- › Available in all upholstery color options
- › Two all metal full extension storage drawers
- › Leveling glides
- › 5-year limited warranty
- › Paper cutter band
- › Paper rod holder

Specifications

Height: 31.5" (80cm)
 Length: 70" (177.8cm)
 Width: 26.5" (67.3cm)
 Storage drawers: 27" W x 16" D x 4.2" H
 (67cm x 41cm x 11cm)
 Net Weight: 132lb (59.9kg)
 Shipping Weight: 185lb (84 kg)
 Shipping Dimensions: 76"L x 30"W x 41"H
 (193cm x 76.2cm x 104.1cm)