

Ultra Comfort Spring Stool with Backrest and Foot Ring

Model 6740

All UMF Medical stools are available in 16 best-in-class seamless upholstery colors with Prefixx® Protective Finish to assist with infection prevention and for full exam room cohesion.



 ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Colors



Key features at a glance

- › Ultra Comfort seamless seat
- › Convenient foot ring
- › Padded backrest (color matched to seat)
- › 5 leg black plastic base
- › 2" (5.1cm) twin wheel casters
- › Air spring height adjustment
- › Available in all UMF Medical upholstery colors
- › 3-year limited warranty

Specifications

Seat diameter: 16" (40.6cm)
 Min. seat height: 16" (40.64cm)
 Max. seat height: 21" (53.34cm)
 Net weight: 22lb (10kg)
 Shipping weight: 27lb (12.2kg)
 Shipping dimensions: L23xW24xH12 (58.4cm x 61cm x 30.5cm)

Accessories

150 – Stool glides – set of 5
 111 – Auto-lock occupied casters