



# Ultra Comfort Spring Stool with Backrest and Foot Ring

## Model 6740

All UMF Medical stools are available in 16 best-in-class seamless upholstery colors with Prefixx® Protective Finish to assist with infection prevention and for full exam room cohesion.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

## Colors



### Key features at a glance

- › Ultra Comfort seamless seat
- › Convenient foot ring
- › Padded backrest (color matched to seat)
- › 5 leg black plastic base
- › 2" (5.1cm) twin wheel casters
- › Air spring height adjustment
- › Available in all UMF Medical upholstery colors
- › 3-year limited warranty

### Specifications

Seat diameter: 16" (40.6cm)  
 Min. seat height: 17.25" (43.8cm)  
 Max. seat height: 22.5" (57.2cm)  
 Net Weight: 22lb (10kg)  
 Shipping Weight: 27lb (12.2kg)  
 Shipping Dimensions: L23xW24xH12 (58.4cm x 61cm x 30.5cm)

### Accessories

150 – Stool glides – set of 5  
 111 – Auto-lock occupied casters



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
 1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com