

Ultra Comfort Air Spring Stool

Model 6749

All UMF Medical stools are available in 16 best-in-class seamless upholstery colors with Prefixx® Protective Finish to assist with infection prevention and for full exam room cohesion.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Colors



Key features at a glance

- › Ultra Comfort seamless seat with anti-microbial guard
- › 5 leg black plastic base
- › 2" (5.1cm) twin wheel casters
- › Air spring height adjustment
- › Available in all UMF upholstery colors
- › 3-year limited warranty

Specifications

Seat diameter: 16" (40.6cm)
 Min. seat height: 17.25" (43.8cm)
 Max. seat height: 22.5" (57.2cm)
 Net Weight: 13lb (5.9kg)
 Shipping Weight: 18lb (8.2kg)
 Shipping Dimensions: L23xW24xH8 (58.4cm x 61cm x 20.3cm)

Accessories

6749F – Fire retardant foam seat cushion
 150 – Stool glides – set of 5
 111 – Auto-lock occupied casters