



Phlebotomy Chair with Foot Operated Pump

Model 8675

The UMF Medical 8675 offers a manual backrest and leg rest—each with its own adjustment lever—for reclining a patient to 10 degrees from horizontal position, as well as a manual foot pedal height adjustment. It has a 360-degree rotation with swivel lock and a 250-lb weight capacity.



Actual colors may vary from the images shown due to monitor settings

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



Key features at a glance

- › 250lb (113.4kg) patient weight capacity
- › Ultra comfort seamless top with antimicrobial guard
- › Sturdy steel frame
- › 360° rotation with swivel lock
- › 10° elevation from horizontal position
- › Manual foot pedal height adjustment
- › Manually adjusted back and leg rest for reclining a patient to a horizontal position
- › Hinged arms
- › 3-year limited warranty

Specifications

Length: 67" (170.2cm) fully extended
 Min. Seat Height: 23" (58.4cm)
 Max. Seat Height: 38" (96.5cm)
 Overall Width: 31.5" (80cm)
 Seat Width: 18" (45.72cm)

Arm Length: 23" (58.4cm)
 Net Weight: 171lb (77.6kg)
 Shipping Weight: 230lb (104.4kg)
 Shipping Dimensions: L44xW34xH64 (111.8cm x 86.4cm x 162.6cm)



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
 1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com