

# Pediatric Examination Table

## Model 5900

UMF Medical Pediatric Exam Table with 2 drawers and a 2-door spacious storage cabinet. Includes digital scale with lbs/oz and kg options, tare weight capability and measuring rule with paddles.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

### Colors



### Key features at a glance

- › Upholstery with Prefixx® Protective Finish
- › Heavy duty roller bearing drawer glides
- › Padded front bumper for infant safety
- › All metal cabinet doors with commercial grade hinges
- › 18" (45.7cm) paper roll holder
- › Digital scale operated by 6 D-cell batteries (not included) with auto shut-off feature
- › Weight options in pounds, ounces and/or kilograms
- › Tare weight capability
- › Scale programmability to KG option only
- › 3-year limited warranty

### Specifications

Overall height: 44.5" (113cm)  
 Work surface height: 35" (88.9cm)  
 Width: 48" (121.9cm)  
 Depth: 25" (63.5cm)  
 Drawers: 18"W x 14" D x 14" H (45cm x 35cm x 10cm)  
 Overall door opening: 40.25" W x 17.25" H (102cm x 44cm)

Net weight: 163lb (74kg)  
 Shipping weight: 233lb (105.7kg)  
 Shipping dimensions: L 60 x W 30 x H 47 (152.4cm x 76.2cm x 119.4cm)  
 Tabletop working dimensions: 48" (121.9cm)

### Accessories

095 – Power cord  
 502 – Double door lock (factory installed only)