

# Bassinet, Two Door with Basket and Mattress

## Model SS8528

UMF Medical's stainless steel bassinets all come standard with the basket (8545), mattress (8546) and swivel casters included. Manufactured on a heavy duty 1" O.D. (2.54cm) tubular stainless steel frame, and basket supports with insert bends, allowing the basket to be moved to horizontal or elevated positions. All bassinets come with a 3-year limited warranty.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › Heavy duty tubular stainless steel framework
- › Full width stainless steel drawer
- › Convenient storage shelf
- › 3" (7.6cm) ball bearing swivel casters
- › Includes basket and mattress pad
- › Name card holder
- › 3-year limited warranty

### Specifications

Height: 35.375" (89.9cm)  
 Width: 31.25" (79.4cm)  
 Depth: 17" (43.2cm)

### Accessories

132 – Swivel casters locking 3" (set of 4)