

24" Base Cabinet with Two Drawers and Two Doors

Model 6014

UMF Medical's 6014 Modular Base Cabinet includes two drawers, two doors and a high-pressure laminate countertop with backsplash.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › High pressure laminate countertop with 4" (10.2cm) backsplash
- › Commercial grade, adjustable, spring assisted door hinges
- › Full extension roller bearing drawer glides
- › Interchangeable drawer and door fronts
- › 3-year limited warranty

Specifications

Base cabinet only (including drawers):	Base cabinet with countertop:
Width: 24" (61cm)	Width: 24" (61cm)
Height: 34.5 (87.6cm)	Height: 35.875" (91.1cm) – does not include height of backsplash
Depth: 18" (45.7cm)	Height: 39.75" (101cm) – includes height of backsplash
	Depth: 19" (48.3cm)

Accessories

- 505 – Single drawer lock (*factory installed*)
- 507 – Additional Keys (x2)